

Recipes:

A. Resist Type: PMMA (+ve) e-beam resist

Supplier order code: Microchem

Resist name: NANO 950 PMMA A4 (4% solid contents)

Molecular weight: 950 K

Solvent: in Anisole

Dehydration bake: 5 minute at 170 deg.

Spinning details:

Speed: 4000rpm, Speed: 3000 rpm (for lift-off)

[Please check the data sheet if one needs some other resist thickness]

Pre bake: 90 sec at 180 deg.

Exposure:

Acc voltage: 10KV

Dose: 80uC/cm²- 200uC/cm² as per your feature size and development time

Write field: 100um

Developer: IPA: MIBK::3:1, developing time: 30 s

Stopper: 2-propanol (IPA) pure, Stopping time: 15 s

Dry: with a low-nitrogen flow

B. Resist Type: (+ve) e-beam resist [Bilayer process]

Supplier order code: Microchem

Resist Name: NANO 950 PMMA A4 (4% solid contents) and CO-polymer EL9

Solvent: in Anisole

Dehydration bake: 5 minute at 170 deg.

Layer1: EL9, Spin @ 3200rpm for 60 Sec

Prebake: 7min @180°C

Layer2: PMMA950K 4%, Spin@2000rpm for 45 sec

Prebake: 2min@180°C

[Please check the data sheet if one needs some other resist thickness]

Exposure:

Acc voltage: 10KV, Apert- 20um

Dose: 100uC/cm²~160uC/cm²

Developer: IPA: MIBK::3:1, developing time: 20s ~90s

Stopper: 2-propanol (IPA) pure, Stopping time: 15s ~ 60s

Dose and development time can be varied depending on the feature sizes

Dry: with a low-nitrogen flow

C. Resist Type: (+ve) e-beam resist [Bilayer process]

Supplier order code: Microchem

Resist Name: NANO 950 PMMA A2 (2% solid contents) and CO-polymer EL9

Solvent: in Anisole

Dehydration bake: 5 minute at 180 deg.

Layer1: EL9, Spin @ 3500rpm for 60 Sec

Prebake: 7min @180°C

Layer2: PMMA950K 2%, Spin@4500rpm for 45 sec

Prebake: 2min@180°C

[Please check the data sheet if one needs some other resist thickness]

Exposure:

Acc voltage: 20KV, Apert-7.5um

Dose: 100uC/cm²~130uC/cm²

Developer: IPA: MIBK::3:1, developing time: 20s

Stopper: 2-propanol (IPA) pure, Stopping time: 15s

Dose and development time can be varied depending on the feature sizes

Dry: with a low-nitrogen flow

D. Resist Type: HSQ (-ve) e-beam resist
Supplier order code: Dow Corning XR -1541 - 6%
Solvent: MIBK
Dehydration bake: 15 minutes at 170 deg.
Speed: 4000rpm

Pre bake: 2 min at 250 deg.

Exposure:

Acc voltage: 10KV

Dose: 200uC/cm²

Write field: 100um

Developer: 25% TMAH , Developing time: 15s

Rinse: Flowing DI-Water, Rinsing time: 30s

Dry: with a low-nitrogen flow.

E. Resist Type: SU8-2000.5 (-ve) e-beam resist
Supplier order code: Microchem SU8-2000.5
Dehydration bake: 5 minutes at 95 deg.
Spinning: Speed: 6000rpm

Pre bake: 65 deg- 1min - ramp to - 95 deg- 1min

Exposure-

Acc voltage : 10KV

Aperture : 7.5um

Dose : 4 uC/cm²

Post Bake: 65 deg- 1min - ramp to- 95 deg- 1min

Developer: SU8 developer- Developing time: 1min

Rinse: IPA- Rinsing time: 1min

Dry: with a low-nitrogen flow.

SU8 Removal

- Soak in PG Remover @ 80 – 30min
- Change container : Again soak in PG remover @ 80 deg - 30min
- Ultra sonication @80deg in PG remover- 5~15 min

Time will depend on the pattern sizes

F. Resist Type: (+ve) e-beam resist

Supplier order code: Zeon Chemicals

Resist Name: ZEP

Solvent: Anisole

Dehydration bake: 2~5 minute at 170 deg.

Layer1: ZEP, Spin @ 2000rpm for 60 Sec

Prebake: 3min @170°C

[Please check the data sheet if one needs some other resist thickness]

Exposure:

Acc voltage: 20KV

Dose: 40uC/cm²

Aperture- 30

Developer: ZED N50, developing time: 60s

Rinse: 2-propanol (IPA) pure, Time: 30s

Post Bake - 100°C – 2 Mins

Dry: with a low-nitrogen flow

G. Resist Type: HSQ (-ve) e-beam resist (Low Temp)

Supplier order code: Dow Corning XR -1541 - 6%

Solvent: MIBK

Dehydration bake: 5 minutes at 170 deg.

Speed: 5000rpm @ 2500 for 60s

Pre bake: 5min at 90 deg.

Exposure:

Acc voltage: 15KV

Dose: 1000 ~ 1500uC/cm²

Write field: 100um

Aperture – 20um

Developer: 25% TMAH , Developing time: 5s

Rinse: Flowing DI-Water, Rinsing time: 30s

Dry: with a low-nitrogen flow.
