

XP-2 High Performance Surface Profiler

Profiler Specifications

System can be used till 100 micron range

Sample Stage Diameter 200mm

Sample Thickness 30mm maximum

Scan Length Range 50mm maximum

X-Y Stage Translation 150mm X 178mm


Vertical Range 400 μ m maximum (200 μ m if up and down step is used)

Vertical Resolution 1 \AA at 10 μ m, 15 \AA at 100 μ m, 62 \AA at 400 μ m

Stylus Tip Radius 2.0 microns

Stylus Force Range .05-10mg (programmable)

Profilometer Operating steps

- Turn on the vacuum pump from the service corridor
- Switch **ON** the main **Power Supply**
- Switch **ON** the computer
- Switch **ON** the power to the **XP Profiler from the back panel**
- **Double-click** on the **XP2 software**  icon located on the desktop
- The main screen is displayed with a dialogue box asking “**Home the System**” before moving on, at the screen prompt, **Left Click** the **Yes** button to continue. The system will move the stage and measurement head stage to their home positions
- A dialogue box appears with a message “**Homing is Complete**”. **Left click** on the **OK** button to continue
- The system is now ready to use
- **Loading the sample:** Open the sample stage cover & place the sample on the stage table and close the sample stage door cover
- **Click** on the Vacuum radio button located on the Main software screen
- Ensure that the stylus is raised enough so as to not hit the sample
- Move the stage in Y direction using the X-Y stage controls or use the Load Position button until the sample is positioned under the measurement head (Stylus)
- Use the X-Y controls to “fine” position the sample for measurement
- The stylus can be raised or lowered at any time by using the mouse and left clicking on the up or down arrows located in the Z stage control panel
- Press the **setup button** located on the Main software screen and edit the parameters according to the requirement
- **Click** on the **Engage** button located on the Z stage control panel

- **Click on Scan** if the stylus is positioned at the required place
- The **Real Time Scan Display** (located on the main screen below the Video window) shows a trace as it is generated

- The Data Display Analysis screen is displayed; the Reference & Measurement cursors appear at their preset position
 - Level Data-** Place the reference (R) and measurement (M) cursors at the data which represents the substrate surface

 - Data width** - the horizontal distance between the R & M cursors along the X axis

 - Delta Height** - the vertical distance between the R & M cursors where they intersect the profile along the Z axis

- For **Saving** data , go to **File→ Save Data File**
- To remove the sample from the instrument, move the Z Stage by using “up” arrow to lift the stylus off the sample surface to a safe distance
- Click on the **LOAD POSITIONS** on the X-Y stage control panel. Open the stage door switch off the vacuum and remove the sample from the stage table
- Close the software and shut down the computer and system