

How to clean and fill your syringe and your bowl?

a) The syringe

- Firstly, **unscrew** the aluminium block syringe.
- Carefully **clean** the syringe and its plunger with ethanol to insure no dust or contamination.
- Then, **clean** it with the liquid to be used (fill it up 3 times).
- **Fill** up your liquid inside the syringe upside down.
- **Insure** no bubble inside the syringe to avoid leak during the measurement. For that, please hold the syringe and **push** the air bubble by pushing the plunger until a drop of liquid comes out of the syringe.
- **Place** your syringe in the aluminium block syringe then screw the two parts of the block. Do **not force**! If a resistance occurs, make sure the small metal rod inside the top block is well rise (unscrew the outside ring if it is not the case).
- **Place** the syringe on the measurement bench, fix it and turn the ring to the top in the direction of clockwise so that the liquid flows.

b) The bowl

- **Clean** it 3 times with ethanol, 3 times with water and 2 times with the liquid to be used.
- **Fill** it with the liquid
- **Place** it under the syringe and make sure the syringe and ensure that it dips into the liquid.