

Usage Policy for PEDOT:PSS

Experimental Overview

- Spin Coating of PEDOT:PSS using LAF (spin coater).
- Annealing of the spin coated layer in N2 ambient on Hot plate.
- Chemical used with PEDOT:PSS
Dimethyl Sulfoxide (DMSO) as a doping chemical in PEDOT:PSS
- Making solution: After adding DMSO in PEDOT:PSS stirring on magnetic stirrer.

Continued...

- Place of experiment: Organic Electronics Lab (LAF spin coater and hot plate)
- Precautionary Measures:
 1. Use safety glass with side shield to avoid eye contact.
 2. Impervious Clothing(Lab coat, Hand gloves, Hair-net)
 3. Proper Exhaust ventilation.
 4. Use mask to avoid ingestion.